

# Stericheck

## Closed Sterility Test System



- Minimizes risk of false positives/negatives
- Meets global regulatory requirements

# Stericheck

## Sterility Test Devices



### Pre-sterilized, Ready to use

- Gamma irradiated
- EO Sterilized

### Validated for

- Sterility
- Microbial Recovery
- Microbial Retention
- Equal Sample Distribution



## Unique Performance Advantage



### Non Coring Needle Design

- To prevent blockage by rubber particles from vial bungs

### Pre-installed Color Coded Clamps

- For easy identification and clamping at both ends of the tube





## Raised Transparent Vent

- To minimize chances of fluid logging and subsequent obstruction of air flow
- For easy visibility in case of accidental fluid logging

## Special 'L' shaped Flow Directors

- To minimize frothing, specially with viscous fluids



## Tough Membrane

- Withstands inadvertent back pressure
- Available with three different polymeric tough membranes

## Archivable Chemical Indicator

- In each Stericheck pack as evidence of EO gas sterilization



## Easy to Peel off, Impervious Blister Pack Cover

- For radiation sterilized Stericheck devices to prevent ingress of  $H_2O_2$  during isolator disinfection

## Customized Needle Connections

- For specific customer needs such as mini vials and cartridges in biopharmaceuticals



# Applications

**Stericheck SVP1**  
Small Volume Parenterals  
Liquid Vials

**Stericheck SVP2**  
Ampules & Collapsible bags

**Stericheck SVP3**  
Small Volume Parenterals  
Vials with soluble powder

**Stericheck SVP4**  
Small Volume Parenterals  
Mini Vials and Cartridges

**Stericheck PC**  
Plastic containers

**Stericheck LVP**  
Large Volume Parenterals  
in glass bottles

**Stericheck PFS**  
Pre-filled Syringes

**SVP Dilutor**  
Antibiotics and difficult to  
dissolve powders

## Ordering Information

### For Stericheck Canisters:

Type				Size		Pore Size		XX	XX	Sterilization		Pack Size	
Type	Code			Dia	Code	Pore Size	Code				Code	Pack Size	Code
	PVDF	Cellulose Nitrate	Nylon-66*	47mm	09	0.45µm	02			EO Sterile	2	10	02
Stericheck SVP1	SV1V	SV1C	SV1N							Gamma Sterile	3		
Stericheck SVP2	SV2V	SV2C	SV2N										
Stericheck SVP3	SV3V	SV3C	SV3N										
Stericheck SVP4	SV4V	SV4C	SV4N										
Stericheck LVP	SLVV	SLVC	SLVN										
Stericheck PC	SPCV	SPCC	SPCN										
Stericheck PFS	SPFV	SPFC	SPFN										

\*Stericheck Devices with Nylon membranes are available as EO sterilized Only

**Example:**

SV1C	09	02	XX	XX	2	02
------	----	----	----	----	---	----

### For SVP Dilutor pack of 20:

EO Sterilized: SVPDXXXXXXXXX209

Gamma Sterilized: SVPDXXXXXXXXX309

### mdi Membrane Technologies INC.

5340 Jaycee Avenue, Suite A, Harrisburg, PA 17112  
Phone: +1-717-412-0943, Fax: +1-717-695-9637  
E-mail: rs@mdimembranetech.com  
Website: www.mdimembranetech.com

### ADVANCED MICRODEVICES PVT. LTD.

20-21, Industrial Area, Ambala Cantt- 133 006, India  
Tel: +91 - 171-2699290, 2699471, Fax: +91-171-2699221  
E-mail: info@mdimembrane.com  
Website: www.mdimembrane.com